## Spectroscopy Lab I

## **ASTR 1010L**

## **Instructor: Larry Krumenaker**



1) Following the instructions given verbally, look at the continuous spectrum light source and make as precise a measurement as you can *for the places where the colors change from one to another*, and *where the red and blue colors change to black*. You may have a partner record YOUR observations but then s/he must look through the spectroscopes and get their own observations. This is not a group observation!

Black-Red	Red-Orange	<b>Orange-Yellow</b>	Yellow-Green	Green-Blue
nm				

	Blue –Indigo	Indigo-Violet	Violet-Black
			nm
2) When you have completed your observations, please enter them in the appropriate columns of the spectroscope spreadsheet on the computer.			

ENTER the LINE NUMBER	on which you put your data in the spreadsheet
here:	(line number)

(There will be a discussion before you do the next part! Please wait!)

3) Go over to Mercury or Helium vapor gas tube and measure the visible lines, recording the *line colors and wavelengths* below. There should be about 8 lines, some strong, some weak. Do not use the Difference column yet.

	Color	Wavelength(nm)	Difference
(Line most towards the violet end of the spectrum)			
(Line most towards the red end of the spectrum)			

- 4) Using the data provided by the instructor, determine which lines you measured and calculate the Difference, as in KNOWN wavelength MINUS YOUR measure. List the data in the Difference column in the table on the other side, being sure to label the numbers with the proper sign (+ or -)
- 5) Graph the values on the grid below. Do NOT connect the points with a line! Notice that we are changing the title from Difference to Correction Factor. You will use the Correction Factor later.



Notes: